

## Anti-Interleukin-6 IL6 Antibody Picoband® PE Conjugated

Catalog Number: PA1700-PE

### About IL6

Interleukin-6 (IL-6) is a protein that in humans is encoded by the IL6 gene. IL-6 is an interleukin that acts as both a pro-inflammatory and anti-inflammatory cytokine. It is secreted by T cells and macrophages to stimulate immune response to trauma, especially burns or other tissue damage leading to inflammation. IL-6 is one of the most important mediators of fever and of the acute phase response. IL-6 is also essential for hybridoma growth and is found in many supplemental cloning media such as briclone. Bowcock et al. (1988) assigned the IL6 gene to chromosome 7p21. By in situ hybridization and Southern blot analysis of mouse-human hybrid cell lines, Sutherland et al. (1988) mapped the IL6 gene to chromosome 7p15.

### Overview

Product Name	Anti-Interleukin-6 IL6 Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P05231

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human IL-6.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 36 Publications Citing This Product

---

1. PubMed ID: 24603891, Sun W, Liu Db, Li Ww, Zhang Li, Long Gx, Wang Jf, Mei Q, Hu Gq. Int J Oncol. 2014 May;44(5):1551-60. Doi: 10.3892/Ijo.2014.2323. Epub 2014 Mar 5. Interleukin-6 Promotes The Migration And Invasion Of Nasopharyngeal Carcinoma Cell Lines And Upregula...
2. PubMed ID: 26990441, Inhibition of Acute Lung Injury by TNFR-Fc through Regulation of an Inflammation-Oxidative Stress Pathway
3. PubMed ID: 26221080, Protective effects of thalidomide on pulmonary injuries in a rat model of paraquat intoxication

Visit [bosterbio.com/anti-il-6-antibody-pa1700-boster.html](http://bosterbio.com/anti-il-6-antibody-pa1700-boster.html) to see all 36 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-6 IL6 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.