

Anti-Interleukin-13 IL13 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PA1701-Fluoro550

About IL13

IL-13, Interleukin 13, is a protein that in humans is encoded by the IL-13 gene. IL-13 is cytokine secreted by many cell types, but especially T helper type 2 (Th2) cells, that is an important mediator of allergic inflammation and disease. The IL-13 gene is mapped to 5q23-q31. IL-13 induces its effects through a multi-subunit receptor that includes the alpha chain of the IL-4 receptor (IL-4Ralpha) and at least one of two known IL-13-specific binding chains. IL-13 acts more prominently as a molecular bridge linking allergic inflammatory cells to the non-immune cells in contact with them, thereby altering physiological function.

Overview

Product Name	Anti-Interleukin-13 IL13 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P35225

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human IL-13.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

7 Publications Citing This Product

1. PubMed ID: 24498391, Li J, Li L, Zuo H, Ke C, Yan B, Wen H, Zhang Y, Yang X. Plos One. 2014 Feb 3;9(2):E87887. Doi: 10.1371/Journal.Pone.0087887. Ecollection 2014. T-Helper Type-2 Contact Hypersensitivity Of Balb/C Mice Aggravated By Dibutyl Phthalate Via Long-Term De...
2. PubMed ID: 26258681, An IL-13 Promoter Polymorphism Associated with Liver Fibrosis in Patients with Schistosoma japonicum
3. PubMed ID: 19739097, SUMF2 interacts with interleukin-13 and inhibits interleukin-13 secretion in bronchial smooth muscle cells

Visit bosterbio.com/anti-il-13-antibody-pa1701-boster.html to see all 7 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-13 IL13 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.