

Anti-Interleukin-13 IL13 Antibody Picoband® PE Conjugated

Catalog Number: PA1701-PE

About IL13

IL-13, Interleukin 13, is a protein that in humans is encoded by the IL-13 gene. IL-13 is cytokine secreted by many cell types, but especially T helper type 2 (Th2) cells, that is an important mediator of allergic inflammation and disease. The IL-13 gene is mapped to 5q23-q31. IL-13 induces its effects through a multi-subunit receptor that includes the alpha chain of the IL-4 receptor (IL-4Ralpha) and at least one of two known IL-13-specific binding chains. IL-13 acts more prominently as a molecular bridge linking allergic inflammatory cells to the non-immune cells in contact with them, thereby altering physiological function.

Overview

Product Name	Anti-Interleukin-13 IL13 Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P35225

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human IL-13.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

7 Publications Citing This Product

1. PubMed ID: 24498391, Li J, Li L, Zuo H, Ke C, Yan B, Wen H, Zhang Y, Yang X. Plos One. 2014 Feb 3;9(2):E87887. Doi: 10.1371/Journal.Pone.0087887. Ecollection 2014. T-Helper Type-2 Contact Hypersensitivity Of Balb/C Mice Aggravated By Dibutyl Phthalate Via Long-Term De...
2. PubMed ID: 26258681, An IL-13 Promoter Polymorphism Associated with Liver Fibrosis in Patients with Schistosoma japonicum
3. PubMed ID: 19739097 , SUMF2 interacts with interleukin-13 and inhibits interleukin-13 secretion in bronchial smooth muscle cells

Visit bosterbio.com/anti-il-13-antibody-pa1701-boster.html to see all 7 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-13 IL13 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.