

Anti-Caspase-8(p18)/CASP8 Antibody Picoband®

Catalog Number: PA1735

About CASP8

Caspase 8 is a caspase protein. It most likely acts upon caspase 3. This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. The human CASP8 gene, whose product is also known as caspase 8 and FLICE, encodes an interleukin-1beta converting enzyme (ICE)-related cysteine protease that is activated by the engagement of several different death receptors. Caspase 8 is immediately recruited to the Fas receptor once it oligomerizes, and its protease activity is crucial for the apoptotic response generated by the resulting death-inducing signaling complex (DISC). This gene contains at least 11 exons spanning approximately 30kb on human chromosome band 2q33-34. This region of human chromosome 2 was previously reported as the location of the CASP10 gene, whose product is closely related to caspase 8. Caspase-8 deficiency in humans is compatible with normal development and shows that caspase-8 has a postnatal role in immune activation of naive lymphocytes.

Overview


Product Name	Anti-Caspase-8(p18)/CASP8 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-Caspase-8 (p18)/CASP8 Antibody catalog # PA1735. Tested in IHC, WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Thimerosal, 0.05mg NaN ₃ . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14790

Technical Details


Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Caspase-8(p18).
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot, and HRP Conjugated anti-Rabbit IgG Super Vision Assay Kit (SV0002-1) for IHC(P).

Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry (Paraffin-embedded Section), 0.5-1ug/ml, Human Western blot, 0.1-0.5ug/ml, Human

Anti-Caspase-8(p18)/CASP8 Antibody Picoband® (PA1735) Images



Anti-Caspase-8(p18) antibody, PA1735, Western blotting
All lanes: Anti Caspase-8(p18) (PA1735) at 0.5ug/ml
Lane 1: A549 Whole Cell Lysate at 40ug
Lane 2: HELA Whole Cell Lysate at 40ug
Predicted bind size: 18KD
Observed bind size: 18KD



Anti-Caspase-8(p18) antibody, PA1735, IHC(P)
IHC(P): Human Tonsil Tissue

16 Publications Citing This Product

1. PubMed ID: 25672487, Zhan W, Hu X, Yi J, An Q, Huang X. Mol Med Rep. 2015 Jun;11(6):4142-8. Doi: 10.3892/Mmr.2015.3326. Epub 2015 Feb 10. Inhibitory Activity Of Apogossypol In Human Prostate Cancer In Vitro And In Vivo.
2. PubMed ID: 28032492, Croton Tiglium Extract Induces Apoptosis via Bax/Bcl-2 Pathways in Human Lung Cancer A549 Cells
3. PubMed ID: 26552967, 2-DG-Regulated RIP and c-FLIP Effect on Liver Cancer Cell Apoptosis Induced by TRAIL

Visit bosterbio.com/anti-caspase-8-p18-antibody-pa1735-boster.html to see all 16 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caspase-8(p18)/CASP8 Antibody

For Research Use Only. Not for use in diagnostic procedures.