

## Anti-Caspase-8(p18)/CASP8 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PA1735-Fluoro550

### About CASP8

Caspase 8 is a caspase protein. It most likely acts upon caspase 3. This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. The human CASP8 gene, whose product is also known as caspase 8 and FLICE, encodes an interleukin-1beta converting enzyme (ICE)-related cysteine protease that is activated by the engagement of several different death receptors. Caspase 8 is immediately recruited to the Fas receptor once it oligomerizes, and its protease activity is crucial for the apoptotic response generated by the resulting death-inducing signaling complex (DISC). This gene contains at least 11 exons spanning approximately 30kb on human chromosome band 2q33-34. This region of human chromosome 2 was previously reported as the location of the CASP10 gene, whose product is closely related to caspase 8. Caspase-8 deficiency in humans is compatible with normal development and shows that caspase-8 has a postnatal role in immune activation of naive lymphocytes.

### Overview

Product Name	Anti-Caspase-8(p18)/CASP8 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q14790

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Caspase-8(p18).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## 16 Publications Citing This Product

1. PubMed ID: 25672487, Zhan W, Hu X, Yi J, An Q, Huang X. Mol Med Rep. 2015 Jun;11(6):4142-8. Doi: 10.3892/Mmr.2015.3326. Epub 2015 Feb 10. Inhibitory Activity Of Apogossypol In Human Prostate Cancer In Vitro And In Vivo.

2. PubMed ID: 28032492, Croton Tiglium Extract Induces Apoptosis via Bax/Bcl-2 Pathways in Human Lung Cancer A549 Cells

3. PubMed ID: 26552967, 2-DG-Regulated RIP and c-FLIP Effect on Liver Cancer Cell Apoptosis Induced by TRAIL

Visit [bosterbio.com/anti-caspase-8-p18-antibody-pa1735-boster.html](http://bosterbio.com/anti-caspase-8-p18-antibody-pa1735-boster.html) to see all 16 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caspase-8(p18)/CASP8 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.