

Anti-Interleukin-4 IL4 Antibody Picoband® HRP Conjugated

Catalog Number: PA1749-HRP

About IL4

Interleukin-4 (IL-4), also known as a B-cell stimulatory factor 1 (BSF1), is an immunomodulatory cytokine, which can inhibit the growth of tumour cells. The human cDNA contains a single open reading frame encoding a protein of 153 amino acids, including a putative signal peptide. IL-4 may act as an autocrine growth factor in pancreatic cancer cells and also give rise to the possibility that cancer-derived IL-4 may suppress cancer-directed immunosurveillance in vivo in addition to its growth-promoting effects, thereby facilitating pancreatic tumor growth and metastasis. The mouse and human genes and their protein products show structural and functional similarities. The human IL-4 gene, which occurs as a single copy in the haploid genome, is mapped on chromosome 5. The standard product used in this kit is recombinant human IL-4, consisting of 130 amino acids with the molecular mass of 14kDa.

Overview

Product Name	Anti-Interleukin-4 IL4 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P05112

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human IL-4.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

7 Publications Citing This Product

1. PubMed ID: 10.1016/j.humpath.2018.04.007, IL-13 is produced by tumor cells in breast implant-associated anaplastic large cell lymphoma: implications for pathogenesis
2. PubMed ID: 23065188, Wang Yz, Feng Xg, Shi Qg, Hao Yl, Yang Y, Zhang Am, Kong Qx. Inflammation. 2013 Apr;36(2):337-45. Doi: 10.1007/S10753-012-9551-5. Silencing Of Mir155 Promotes The Production Of Inflammatory Mediators In Guillain-Barr?? Syndrome In Vitro.
3. PubMed ID: 23011634, Tan Y, Liu H, Yang H, Wang L, Qin X. J Physiol Sci. 2013 Jan;63(1):63-9. Doi: 10.1007/S12576-012-0231-X. Epub 2012 Sep 26. An Inactivated Pseudomonas Aeruginosa Medicament Inhibits Airway Allergic Inflammation And Improves Epithelial Functions.

Visit bosterbio.com/anti-il-4-antibody-pa1749-boster.html to see all 7 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-4 IL4 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.