

Anti-RIP2/RIPK2 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA1861-Fluoro488

About RIPK2

RIPK2 (Receptor-interacting serine/threonine-protein kinase 2), also known as CARD3, CARDIAK, RICK, RIP2, is an enzyme that in humans is encoded by the RIPK2 gene. It has 540-amino acid protein in length. Northern blot analysis revealed that RICK is expressed in various human tissues as 2.5- and 1.8-kb mRNAs that differ due to alternative polyadenylation. RICK is a novel kinase that may regulate apoptosis induced by the FAS receptor pathway. This gene encodes a member of the receptor-interacting protein (RIP) family of serine/threonine protein kinases. The encoded protein contains a C-terminal caspase recruitment domain (CARD), and is a component of signaling complexes in both the innate and adaptive immune pathways. It is a potent activator of NF-kappa B and inducer of apoptosis in response to various stimuli, CARDIAK (CARD-containing ICE-associated kinase) specifically interacted with the CARD of ICE/caspase-1, and this interaction correlated with the processing of pro-caspase-1 and the formation of the active caspase-1 p20 subunit.

Overview

Product Name	Anti-RIP2/RIPK2 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O43353

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human RIP2, different from the related rat and mouse sequences by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm

	Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RIP2/RIPK2 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.