

Anti-Glutaredoxin 2/GLRX2 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PA1885-Fluoro550

About GLRX2

GLRX2 (Glutaredoxin-2) also known as Glutaredoxin-2, mitochondrial or GRX2, is a protein that in humans is encoded by the GLRX2 gene. Glutaredoxins (e.g., GLRX) are a family of glutathione-dependent hydrogen donors that participate in a variety of cellular redox reactions. By sequence analysis, Lundberg et al. (2001) and Gladyshev et al. (2001) identified the GLRX2 gene within a clone mapping to chromosome 1q31.2-q31.3. Lundberg et al. (2001) determined that the GLRX2 gene contains 5 exons spanning about 9.6 kb. The GLRX2B transcript uses a first exon (exon 1B) located upstream from the first exon used by the GLRX2A transcript (exon 1A), suggesting that alternative splicing generates the isoforms. Lundberg et al. (2001) assayed reductase activity in recombinant proteins of both GLRX2 isoforms. They found that both have GSH-dependent dehydroascorbate reductase activity and 2-hydroxyethyl disulfide reductase activity.

Overview

Product Name	Anti-Glutaredoxin 2/GLRX2 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NS18

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Glutaredoxin 2, different from the related rat sequence by two amino acids, and from the related mouse sequence by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550

	Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 10.1016/j.envint.2020.105654, Enhanced expression of nicotinamide nucleotide transhydrogenase (NNT) and its role in a human T cell line continuously exposed to asbestos

Visit bosterbio.com/anti-glutaredoxin-2-antibody-pa1885-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glutaredoxin 2/GLRX2 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.