

## Anti-Neurosecretory protein VGF VGF Antibody Picoband® PE Conjugated

Catalog Number: PA1984-PE

### About VGF

VGF (VGF, Nerve Growth Factor-Inducible), is a protein and neuropeptide that may play a role in regulating energy homeostasis, metabolism. Canu et al. (1997) noted that rat Vgf encodes a predicted 70-kD polypeptide that shares similarities with the secretogranin/chromogranin family and is found in the secretory granules of subsets of neurons and endocrine cells. By fluorescence in situ hybridization, Canu et al. (1997) assigned the VGF gene to 7q22. Canu et al. (1997) demonstrated that the single-copy human VGF gene spans 6 kb of genomic DNA and contains 2 exons. The entire VGF protein is encoded by exon 2, while exon 1 contains only 5-prime untranslated sequence. The structural organization of the human gene is similar to that described for the rat VGF gene (Salton et al., 1991), and both the translated and the untranslated regions show a high degree of sequence homology to the rat gene.

### Overview

Product Name	Anti-Neurosecretory protein VGF VGF Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O15240

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human VGF, different from the related rat and mouse sequences by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE

	Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Neurosecretory protein VGF VGF Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.