

## Anti-Protein Wnt-2b Wnt2b Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA1985-Fluoro488

### About WNT2B

WNT2B (Wingless-Type MMTV Integration Site Family Member 2B), is a protein that in humans is encoded by the WNT2B gene. This gene encodes a member of the wingless-type MMTV integration site (WNT) family of highly conserved, secreted signaling factors. Katoh et al. (1996) used fluorescence in situ hybridization to map the WNT13 gene to chromosome 1p13. WNT family members function in a variety of developmental processes including regulation of cell growth and differentiation and are characterized by a WNT-core domain. This gene may play a role in human development as well as human carcinogenesis.

### Overview

Product Name	Anti-Protein Wnt-2b Wnt2b Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q93097

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Wnt2b, identical to the related rat and mouse sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

---

1. PubMed ID: 33359712, An R,Wang X,Yang L,Zhang J,Wang N,Xu F,Hou Y,Zhang H,Zhang L. Polystyrene microplastics cause granulosa cells apoptosis and fibrosis in ovary through oxidative stress in rats.Toxicology.2020 Dec 24;449:152665. doi:10.1016/j.tox.2020.152665.Epub ahead of p

2. PubMed ID: 32463570, Xu T,Pan L,Li L,Hu S,Zhou H,Yang C,Yang J,Li H,Liu Y,Meng X,Li J.MicroRNA-708 modulates Hepatic Stellate Cells activation and enhances extracellular matrix accumulation via direct targeting TMEM88. J Cell Mol Med. 2020 Jul;24(13):7127-7140.doi: 10.1111/jc

Visit [bosterbio.com/anti-wnt2b-antibody-pa1985-boster.html](https://bosterbio.com/anti-wnt2b-antibody-pa1985-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein Wnt-2b Wnt2b Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.