

## Anti-TORC1/CRTC1 Antibody Picoband® HRP Conjugated

Catalog Number: PA1987-HRP

### About CRTC1

CRTC1 (CREB-Regulated Transcription Coactivator 1), also known as MECT1, TORC1 or KIAA0616, is a protein that in humans is encoded by the CRTC1 gene. By sequence analysis, Tonon et al. (2003) mapped the CRTC1 gene to chromosome 19p13. Kovacs et al. (2007) found that Torc1 was expressed in adult mouse brain and cultured neurons, and that it translocated to the nucleus upon concomitant activation of calcium and cAMP signaling pathways. Mair et al. (2011) demonstrated that CRTC1 is a direct AMPK target, and interacts with the CREB homolog-1 (CRH1) transcription factor in vivo. The prolongevity effects of activating AMPK or deactivating calcineurin decrease CRTC1 and CRH1 activity and induce transcriptional responses similar to those of CRH1-null worms.

### Overview

Product Name	Anti-TORC1/CRTC1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6UUV9

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human TORC1, different from the related rat and mouse sequences by two amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users.

Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.  
ELISA, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TORC1/CRTC1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.