

## Anti-Nucleoplasmin-2 NPM2 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PA1991-Fluoro550

### About NPM2

NPM2 (Nucleophosmin/Nucleoplasmin Family Member 2), Burns et al. (2003) mapped the mouse Npm2 gene to chromosome 14. Burns et al. (2003) found that female Npm2-null mice showed defects in preimplantation embryo development, with abnormalities in oocyte and early embryonic nuclei.

### Overview

Product Name	Anti-Nucleoplasmin-2 NPM2 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q86SE8

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human NPM2.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nucleoplasmin-2 NPM2 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.