

## Anti-Flotillin 1/FLOT1 Antibody Picoband® FITC Conjugated

Catalog Number: PA2033-FITC

### About FLOT1

FLOT1 (Flotillin 1), is a protein that in humans is encoded by the FLOT1 gene. The International Radiation Hybrid Mapping Consortium mapped the FLOT1 gene to chromosome 6. Bickel et al. (1997) found that mouse Flot1 behaves as a resident integral membrane protein of caveolae. It consistently copurified with Flot2 and with caveolin-1 in the purification of caveolin-rich membranes. Hazarika et al. (1999) found that stable transfection of Flot1, which they called ESA/flotillin-2, in COS-1 cells induced filopodia formation and changed the epithelial morphology to that of neuronal cells. Santamaria et al. (2005) found that prostate tumor overexpressed gene-1 interacted with flotillin-1 in detergent-insoluble membrane fractions. Flotillin-1 colocalized with PTOV1 at the plasma membrane and in the nucleus, and it entered the nucleus concomitant with PTOV1 shortly before initiation of S phase.

### Overview

Product Name	Anti-Flotillin 1/FLOT1 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O75955

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Flotillin 1, different from the related rat and mouse sequences by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 28522298, Shi, X., Li, W., Liu, H., Yin, D., & Zhao, J. (2017). beta-cyclodextrin induces the differentiation of resident cardiac stem cells to cardiomyocytes through autophagy. *Biochimica et Biophysica Acta (BBA) - Molecular Cell Research*, 1864(8), 1425-1434....

Visit [bosterbio.com/anti-flotillin-1-antibody-pa2033-boster.html](http://bosterbio.com/anti-flotillin-1-antibody-pa2033-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Flotillin 1/FLOT1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.