

Anti-FMO5 Antibody Picoband®

Catalog Number: PA2108

About FMO5

Dimethylaniline monooxygenase[N-oxide-forming] 5 also known as FMO 5 is an enzyme that in humans is encoded by the FMO5 gene. By fluorescence in situ hybridization, this gene was assigned to human chromosome 1q21.1. FMO5 is not an efficient drug-metabolizing enzyme and that it may have an alternative physiologic role. FMO5 transcripts play an important role in progesterone-regulated in breast cancer cells specifically under the control of the progesterone receptor B-isoform.

Overview


Product Name	Anti-FMO5 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-FMO5 Antibody catalog # PA2108. Tested in WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Thimerosal, 0.05mg NaN ₃ . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49326

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human FMO5, different from the related rat sequence by one amino acid, and different from the related mouse sequence by two amino acids.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human

Anti-FMO5 Antibody Picoband® (PA2108) Images



Western blot analysis of FMO5 using anti-FMO5 antibody (PA2108). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human HepG2 whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-FMO5 antigen affinity purified polyclonal antibody (Catalog # PA2108) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for FMO5 at approximately 60 kDa. The expected band size for FMO5 is at 60 kDa.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FMO5 Antibody

For Research Use Only. Not for use in diagnostic procedures.