

Anti-CDC5L Antibody Picoband®

Catalog Number: PA2123

About CDC5L

Cell division cycle 5-like protein, also called PCDC5RP or CDC5-like, is a protein that in humans is encoded by the CDC5L gene. The CDC5L gene was mapped to chromosome 6p21 by fluorescence in situ hybridization. This DNA-binding protein involved in cell cycle control. It may act as a transcription activator. The Component of the PRP19-CDC5L complex forms an integral part of the spliceosome and is required for activating pre-mRNA splicing. The role of CDC5 (Polo-like kinase) in regulating Rho1 is likely to be relevant to cytokinesis and asymmetric cell division in other organisms.

Overview


Product Name	Anti-CDC5L Antibody Picoband®
Reactive Species	Human, Rat
Description	Boster Bio Anti-CDC5L Antibody catalog # PA2123. Tested in WB applications. This antibody reacts with Human, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Thimerosal, 0.05mg NaN ₃ . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99459

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human CDC5L, different from the related rat and mouse sequences by four amino acids.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human, Rat

Anti-CDC5L Antibody Picoband® (PA2123) Images



Western blot analysis of CDC5L using anti-CDC5L antibody (PA2123). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 50ug of sample under reducing conditions. Lane 1: human HEK293 whole cell lysates Lane 2: human THP-1 whole cell lysates Lane 3: rat C6 whole cell lysates After Electrophoresis, proteins were transferred to a Nitrocellulose membrane at 150mA for 50-90 minutes. Blocked the membrane with 5% Non-fat Milk/ TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-CDC5L antigen affinity purified polyclonal antibody (Catalog # PA2123) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:10000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for CDC5L at approximately 100KD. The expected band size for CDC5L is at 92KD.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CDC5L Antibody

For Research Use Only. Not for use in diagnostic procedures.