

## Anti-P2X purinoceptor 7 P2RX7 Antibody Picoband®

Catalog Number: PA2129

### About P2rx7

P2X purinoceptor 7, also called P2X7, is a protein that in humans is encoded by the P2RX7 gene. The product of this gene belongs to the family of purinoceptors for ATP. P2RX7 was mapped to chromosome 12q24 by fluorescence in situ hybridization. It is a receptor for ATP that acts as a ligand-gated ion channel. This gene is responsible for ATP-dependent lysis of macrophages through the formation of membrane pores permeable to large molecules. It could function in both fast synaptic transmission and the ATP-mediated lysis of antigen-presenting cells.

### Overview


Product Name	Anti-P2X purinoceptor 7 P2RX7 Antibody Picoband®
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-P2X purinoceptor 7 P2RX7 Antibody catalog # PA2129. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg Thimerosal, 0.05mg NaN3. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q64663

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of rat P2RX7, different from the related mouse sequence by one amino acid, and from the related human sequence by two amino acids.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human, Rat, Mouse

## Anti-P2X purinoceptor 7 P2RX7 Antibody Picoband® (PA2129) Images



Anti-P2RX7 antibody, PA2129, Western blottingWB: U87 Cell  
Lysate

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-P2X purinoceptor 7 P2RX7 Antibody

For Research Use Only. Not for use in diagnostic procedures.