

Anti-Olig1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PA2213-Fluoro594

About OLIG1

Oligodendrocyte transcription factor 1 (OLIG1), also known as class B basic helix-loop-helix protein 6 (BHLHB6) is a protein that in humans is encoded by the OLIG1 gene. This gene is mapped to 21q22.11. This gene is a member of the oligodendrocyte lineage gene family which encodes basic helix-loop-helix transcription factors. Studies in mice show that it is coexpressed with Olig2 in neural progenitors, and has a role in the development and maturation of oligodendrocytes. OLIG1 gene promotes formation and maturation of oligodendrocytes, especially within the brain. It cooperates with OLIG2 to establish the pMN domain of the embryonic neural tube.

Overview

Product Name	Anti-Olig1 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8TAK6

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Olig1, different from the related mouse and rat sequences by four amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Olig1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.