

Anti-HNF6/ONECUT1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA2214-Fluoro488

About ONECUT1

One cut homeobox 1 (ONECUT1), also called Hepatocyte nuclear factor 6 (HNF6) is found strong expression in liver and lower expression in testis and skin. ONECUT1 gene is mapped to 15q21.3. This gene encodes a member of the Cut homeobox family of transcription factors. Expression of the encoded protein is enriched in the liver, where it stimulates transcription of liver-expressed genes, and antagonizes glucocorticoid-stimulated gene transcription. This gene may influence a variety of cellular processes including glucose metabolism, cell cycle regulation, and it may also be associated with cancer.

Overview

Product Name	Anti-HNF6/ONECUT1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9UBC0

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human HNF6, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HNF6/ONECUT1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.