

## Anti-Diubiquitin/UBD Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA2222-Fluoro488

### About UBD

Ubiquitin D, also called UBD or FAT10 is a protein that in humans is encoded by the UBD gene. This gene is mapped to 6p22.1. Ubiquitin-like protein modifier can be covalently attached to target protein and subsequently leads to their degradation by the 26S proteasome in a NUB1L-dependent manner. This gene may be involved in dendritic cell (DC) maturation, the process by which immature dendritic cells differentiate into fully competent antigen-presenting cells that initiate T-cell responses. It may be involved in the formation of aggresomes when proteasome is saturated or impaired. This gene mediates apoptosis in a caspase-dependent manner, especially in renal epithelium and tubular cells during renal diseases such as polycystic kidney disease and Human immunodeficiency virus (HIV)-associated nephropathy (HIVAN). It probably functions as a survival factor.

### Overview

Product Name	Anti-Diubiquitin/UBD Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O15205

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Diubiquitin, different from the related mouse and rat sequences by five amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Diubiquitin/UBD Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.