

Anti-CAT1/SLC7A1 Antibody Picoband® Biotin Conjugated

Catalog Number: PA2245-Biotin

About SLC7A1

High affinity cationic amino acid transporter 1, also called solute carrier family 7 (cationic amino acid transporter, γ^+ system), member 1 (SLC7A1) or ATRC1 is a protein that in humans is encoded by the SLC7A1 gene. This gene is mapped to 13q12.3. It acts as a cationic amino acid transporter and as an ecotropic retrovirus receptor and plays a role in the Na^+ -independent γ^+ cationic amino acid transport system. This gene is high-affinity, low capacity permease involved in the transport of the cationic amino acids (arginine, lysine and ornithine) in non-hepatic tissues.

Overview

Product Name	Anti-CAT1/SLC7A1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na_2HPO_4 , 0.02% NaN_3 .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P30825

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human CAT1, different from the related mouse and rat sequences by four amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CAT1/SLC7A1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.