

## Anti-Vasopressin V1a receptor AVPR1A Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PA2248-Fluoro647

### About Avpr1a

Arginine vasopressin receptor 1A (officially called AVPR1A) is one of the three major receptor types for arginine vasopressin (AVPR1B and AVPR2 being the others). It belongs to the subfamily of G-protein coupled receptors which includes AVPR1B, V2R and OXT receptors. This gene is mapped to 12q14.2. AVPR1A is present throughout the brain, as well as in the periphery, in the liver, kidney, and vasculature. The protein encoded by this gene acts as receptor for arginine vasopressin. Its activity is mediated by G proteins which stimulate a phosphatidylinositol-calcium second messenger system. The receptor mediates cell contraction and proliferation, platelet aggregation, release of coagulation factor and glycogenolysis.

### Overview

Product Name	Anti-Vasopressin V1a receptor AVPR1A Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P30560

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of rat AVPR1A, different from the related human and mouse sequences by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm

	Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Vasopressin V1a receptor AVPR1A Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.