

Anti-TNFAIP8L3 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PA2257-Fluoro647

About TNFAIP8L3

TNFAIP8L3 belongs to the TNFAIP8 family. It acts as a negative mediator of apoptosis and may play a role in tumor progression. This gene suppresses the TNF-mediated apoptosis by inhibiting caspase-8 activity but not the processing of procaspase-8, subsequently resulting in inhibition of BID cleavage and caspase-3 activation.

Overview

Product Name	Anti-TNFAIP8L3 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5GJ75

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human TNFAIP8L3, different from the related mouse and rat sequences by two amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFAIP8L3 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.