

Anti-APG5L/ATG5 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA2260-Fluoro488

About ATG5

Autophagy protein 5 is a protein that in humans is encoded by the ATG5 gene. It is also known as APG5 or ASP, and this gene is mapped to 6q21. It is found that knockdown of ATG5 in hepatocytes increased triglyceride levels with oleate or a second endogenous stimulus for triglyceride formation. These hepatocytes with ATG5 knockdown also had increased lipid droplet number and size. ATG5 is an E3 ubiquitin ligase which is necessary for autophagy due to its role in autophagosome elongation. It is activated by ATG7 and forms a complex with ATG12 and ATG16L1. This complex is necessary for LC3-1 conjugation to PE to form LC3-II.

Overview

Product Name	Anti-APG5L/ATG5 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9H1Y0

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human APG5L, identical to the related mouse sequence, and different from the related rat sequence by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

7 Publications Citing This Product

1. PubMed ID: -, Shen W, Jia N, Miao J, Chen S, Zhou S, Meng P, Zhou X, Tang L, Zhou L: Penicillium B Protects against Cisplatin-Induced Renal Tubular Cell Apoptosis through Activation of AMPK-Induced Autophagy and Mitochondrial Biogenesis. *Kidney Dis* 2021. doi:10.1159/000514657

2. PubMed ID: 29113167, Yang J, Sheng S, Yang Q, Li L, Qin S, Yu S, Zhang X. *Oncol Lett.* 2017 Nov;14(5):5333-5339. doi: 10.3892/ol.2017.6857. Epub 2017 Aug 31. Endocan silencing induces programmed cell death in hepatocarcinoma

3. PubMed ID: 26459718, Anti-autophagic and anti-apoptotic effects of memantine in a SH-SY5Y cell model of Alzheimer's disease via mammalian target of rapamycin-dependent and -independent pathways

Visit bosterbio.com/anti-apg5l-atg5-antibody-pa2260-boster.html to see all 7 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-APG5L/ATG5 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.