

Anti-MEK4/MAP2K4 Antibody Picoband® Cy3 Conjugated

Catalog Number: PA2275-Cy3

About MAP2K4

Dual specificity mitogen-activated protein kinase kinase 4 (MAP2K4), also called SEK1 or JNKK1, is an enzyme that in humans is encoded by the MAP2K4 gene. It is mapped to 17p12. This gene encodes a dual specificity protein kinase that belongs to the Ser/Thr protein kinase family. This kinase is a direct activator of MAP kinases in response to various environmental stresses or mitogenic stimuli. It has been shown to activate MAPK8/JNK1, MAPK9/JNK2, and MAPK14/p38, but not MAPK1/ERK2 or MAPK3/ERK1. This kinase is phosphorylated, and thus activated by MAP3K1/MEKK. MAP2K4 was a specific activator of JNK1, JNK2, and p38, and it has been shown to interact with FLNC, MAPK8, MAPK8IP3 and AKT1.

Overview

| | |
|----------------------|--|
| Product Name | Anti-MEK4/MAP2K4 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P45985 |

Technical Details

| | |
|------------------|--|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human MEK4, different from the related mouse and rat sequences by one amino acid. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MEK4/MAP2K4 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.