

Anti-SOX7 Antibody Picoband® APC Conjugated

Catalog Number: PA2285-APC

About SOX7

SOX7 belongs to SOX gene family and SOX proteins are transcription factors with critical roles in the regulation of diverse developmental processes. This gene is mapped to 8p23.1. SOX7 gene contains at least 2 exons. In cotransfected 293 cells, Sox7 reduced Wnt (see WNT1)/beta-catenin (see CTNNB1)-stimulated transcription. SOX7 is a potent activator of FGF3 transcription. It not only plays a role in the transcriptional regulation of genes expressed in the hemogenic endothelium but also blocks further differentiation into blood precursors. And it may be required for the survival of both hematopoietic and endothelial precursors during specification.

Overview

Product Name	Anti-SOX7 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BT81

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human SOX7, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SOX7 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.