

Anti-TMEM16A/ANO1 Antibody Picoband® APC Conjugated

Catalog Number: PA2290-APC

About ANO1

Anoctamin-1 (ANO1), also known as oral cancer overexpressed 2 (ORAOV2) or tumor-amplified and overexpressed sequence 2 (TMEM16A), is a protein that in humans is encoded by the ANO1 gene. This gene belongs to a family of membrane proteins containing 8 transmembrane segments, and it is mapped to 11q13.3. ANO1 is a candidate calcium-activated chloride channel that mediates receptor-activated chloride currents in diverse physiologic processes, and it is thought to be responsible for a voltage-sensitive calcium-activated chloride current. Its overexpression was reported in esophageal squamous cell carcinoma and breast cancer progression Crofelemer, an antidiarrhoeal, inhibits this channel. ANO1 has eight transmembrane domains, its pore is large and non-selective, allowing other negatively charged species to permeate.

Overview

Product Name	Anti-TMEM16A/ANO1 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5XXA6

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human TMEM16A, different from the related mouse and rat sequences by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 10.1096/fasebj.2020.34.s1.06115, Intestinal TMEM16A function as a luminal chloride channel

Visit bosterbio.com/anti-tmem16a-antibody-pa2290-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TMEM16A/ANO1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.