

## Anti-CaV1.3/CACNA1D Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA2293-Fluoro488

### About CACNA1D

Cav1.3, also known as the calcium channel, voltage-dependent, L type, alpha 1D subunit (CACNA1D), is a human gene. It is mapped to 3p21.1. Voltage-dependent calcium channels mediate the entry of calcium ions into excitable cells, and are also involved in a variety of calcium-dependent processes, including muscle contraction, hormone or neurotransmitter release, and gene expression. Alpha-1D subunit can mediate DHP-sensitive, high voltage-activated, long-lasting calcium channel activity. CACNA1D can form L-type calcium channels with negative activation thresholds which is essential for normal auditory function and controlling of cardiac pacemaker activity.

### Overview

Product Name	Anti-CaV1.3/CACNA1D Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q01668

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human CaV1.3, different from the related mouse and rat sequences by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

---

1. PubMed ID: 28656274, Guo, Y., Cui, L., Jiang, S., Zhang, A., & Jiang, S. (2017). Proteomics of acute heart failure in a rat post-myocardial infarction model. *Molecular Medicine Reports*, 16(2), 1946-1956. Advance online publication. doi: 10.3892/mmr.2017.6820

Visit [bosterbio.com/anti-cav1-3-antibody-pa2293-boster.html](http://bosterbio.com/anti-cav1-3-antibody-pa2293-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CaV1.3/CACNA1D Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.