

Anti-Kv1.1 potassium channel/KCNA1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PA2296-Fluoro594

About KCNA1

Potassium voltage-gated channel subfamily A member 1, also known as Kv1.1, is a shaker related voltage-gated potassium channel that in humans is encoded by the KCNA1 gene. It is mapped to 12p13.32. The protein functions as a potassium selective channel through which the potassium ion may pass through in consensus with the electrochemical gradient. The N-terminus of the channel is associated with beta subunits that can modify the inactivation properties of the channel as well as affect expression levels. The C-terminus of the channel is complexed to a PDZ domain protein that is responsible for channel targeting.

Overview

Product Name	Anti-Kv1.1 potassium channel/KCNA1 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q09470

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Kv1.1 potassium channel, different from the related mouse sequence by two amino acids, and from the related rat sequence by four amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 10.1002/jnr.25001, Adult-onset hypothyroidism causes mechanical hypersensitivity due to peripheral nerve hyperexcitability based on voltage-gated potassium channel downregulation in male mice

Visit bosterbio.com/anti-kv1-1-potassium-channel-antibody-pa2296-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Kv1.1 potassium channel/KCNA1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.