

Anti-TRPC4 Antibody Picoband® APC Conjugated

Catalog Number: PA2307-APC

About TRPC4

The short transient receptor potential channel 4 (TRPC4), also known as Trp-related protein 4, is a protein that in humans is encoded by the TRPC4 gene. TRPC4 is a member of the transient receptor potential cation channels. It is mapped to 13q13.3. This protein forms a non-selective calcium-permeable cation channel that is activated by Gq-coupled receptors and tyrosine kinases, and plays a role in multiple processes including endothelial permeability, vasodilation, neurotransmitter release and cell proliferation. TRPC4 also has an essential role in endothelial-dependent regulation of vascular tone, endothelial permeability, and neurotransmitter release from thalamic interneurons.

Overview

Product Name	Anti-TRPC4 Antibody Picoband® APC Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9UBN4

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human TRPC4.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRPC4 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.