

## Anti-ADAMTS13 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB10028-Fluoro647

### About ADAMTS13

ADAMTS13 (a disintegrin and metalloproteinase with athrombospondin type 1 motif, member 13), also known as VWFCP, is a zinc-containing metalloprotease enzyme that cleaves von Willebrand factor (vWf), a large protein involved in blood clotting. It is secreted in blood and degrades large vWf multimers, decreasing their activity. This gene encodes a member of a family of proteins containing several distinct regions, including a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. This gene is mapped to 9q34.

### Overview

Product Name	Anti-ADAMTS13 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q76LX8

### Technical Details

Immunogen	E. coli-derived human ADAMTS13 recombinant protein (Position: A299-R488). Human ADAMTS13 shares 85.2% amino acid (aa) sequence identity with mouse ADAMTS13.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADAMTS13 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.