

## Anti-AKR1B1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB10035-Fluoro488

### About AKR1B1

Aldo-keto reductase family 1, member B1 (aldose reductase), also known as AR, is an enzyme that in humans is encoded by the AKR1B1 gene. This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. This member catalyzes the reduction of a number of aldehydes, including the aldehyde form of glucose, and is thereby implicated in the development of diabetic complications by catalyzing the reduction of glucose to sorbitol.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-AKR1B1 Antibody Picoband® Fluoro488 Conjugated  |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | P15121   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E. coli-derived human AKR1B1 recombinant protein (Position: L228-F316). Human AKR1B1 shares 87.5% amino acid (aa) sequence identity with both mouse and rat AKR1B1. |
| Cross Reactivity    | No cross-reactivity with other proteins   |
| Isotype             | Rabbit IgG  |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | Fluoro488<br>Excitation Wavelength: 488 nm<br>Emission Wavelength: 515-545 nm   |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.  |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AKR1B1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.