

## Anti-AQP11 Antibody Picoband® APC Conjugated

Catalog Number: PB10044-APC

### About AQP11

AQP11 has a unique asparagine-proline-alanine (NPA) box distinct from those of other AQPs, suggesting a different pore structure and function. Using Northern blot analysis, they detected highest expression of mouse Aqp11 in testis, followed by liver and kidney. Expression was much weaker in heart, brain, and muscle. Western blot analysis of mouse kidney membrane fractions detected Aqp11 at an apparent molecular mass of 26 kD, lower than the calculated molecular mass of 30 kD. Immunohistochemical analysis localized Aqp11 to mouse renal proximal tubule cells, where it showed a perinuclear distribution. Fluorescence-tagged Aqp11 localized with an endoplasmic reticulum marker.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-AQP11 Antibody Picoband® APC Conjugated   |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | Q8NBQ7   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthetic peptide corresponding to a sequence at the N-terminus of human AQP11, different from the related mouse sequence by two amino acids, and from the related rat sequence by three amino acids. |
| Cross Reactivity    | No cross-reactivity with other proteins.  |
| Isotype             | Rabbit IgG  |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | APC<br>Excitation Wavelength: 633-647 nm<br>Emission Wavelength: 660 nm   |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.  |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AQP11 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.