

## Anti-Ki67/MKI67 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9026-Fluoro488

### About MKI67

Ki-67 (Proliferation-related Ki-67 antigen), also known as MKI67 or KIA, is a protein that in humans is encoded by the MKI67 gene. From study of a panel of human-rodent somatic cell hybrids, it has been demonstrated that a gene involved in the expression of the MKI67 antigen is located on chromosome 10. By in situ hybridization, Fonatsch et al. (1991) regionalized the MKI67 gene to chromosome 10q25-qter. By FISH, Traut et al. (1998) mapped the mouse Mki67 gene to chromosome 7F3-F5. Antigen KI-67 is a nuclear protein that is associated with and may be necessary for cellular proliferation. Furthermore it is associated with ribosomal RNA transcription. Inactivation of antigen KI-67 leads to inhibition of ribosomal RNA synthesis.

### Overview

Product Name	Anti-Ki67/MKI67 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P46013

### Technical Details

Immunogen	E.coli-derived human Ki67 recombinant protein (Position: K2860-I3256).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 62 Publications Citing This Product

1. PubMed ID: PMID:31966612, Clinical biocharacterization of immunophenotype in hepatocellular carcinoma patients
2. PubMed ID: 10.1186/s12935-021-01775-5, Regulator of cullins-1 (ROC1) negatively regulates the Gli2 regulator SUFU to activate the hedgehog pathway in bladder cancer
3. PubMed ID: 10.2485/jhtb.30.273, Elevated CREPT Expression Enhances the Progression of Salivary Gland Adenoid Cystic Carcinoma

Visit [bosterbio.com/anti-ki67-picoband-trade-antibody-pb9026-boster.html](https://bosterbio.com/anti-ki67-picoband-trade-antibody-pb9026-boster.html) to see all 62 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ki67/MKI67 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.