

## Anti-CNTF Antibody Picoband® Cy3 Conjugated

Catalog Number: PB9032-Cy3

### About Cntf

Ciliary neurotrophic factor (CNTF) is a potent polypeptide hormone whose actions appear to be restricted to the nervous system where it promotes survival, neurotransmitter synthesis and neurite outgrowth in certain neuronal populations. The mouse CNTF gene is on mouse chromosome 19 and that its expression is unaffected in the mouse neurologic mutant wobbler, a form of spinal muscular atrophy. The CNTF protein is highly conserved in evolution. The protein is a potent survival factor for neurons and oligodendrocytes, and it may be involved in reducing tissue destruction during inflammatory attacks. CNTF is thought to act centrally by inducing hypothalamic neurogenesis to modulate food intake and peripherally by altering hepatic gene expression, in a manner similar to that of leptin.

### Overview

Product Name	Anti-CNTF Antibody Picoband® Cy3 Conjugated
Reactive Species	Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P51642

### Technical Details

Immunogen	E.coli-derived mouse CNTF recombinant protein (Position: A2-M198). Mouse CNTF shares 83% and 95% amino acid (aa) sequences identity with human and rat CNTF, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CNTF Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.