

## Anti-LIF Antibody Picoband® Biotin Conjugated

Catalog Number: PB9036-Biotin

### About LIF

LIF is a pleiotropic cytokine produced at the maternal-fetal interface which has been shown to play an essential role in implantation in mice. This gene is mapped to 22q11-q12.2, between the Philadelphia translocation BCR gene and the breakpoint of the translocation in cell line GM2324 at 22q12.2. LIF is produced in high amounts by the human endometrium and the trophoblast itself, and LIF receptors are present on cytotrophoblast cells. LIF could, thus, play a role in modulating HLA-G production and immune tolerance at the maternal-fetal interface.

### Overview

Product Name	Anti-LIF Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P15018

### Technical Details

Immunogen	E.coli-derived human LIF recombinant protein (Position: S23-F202). Human LIF shares 78% and 82% amino acid (aa) sequences identity with mouse and rat LIF, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

## 10 Publications Citing This Product

1. PubMed ID: 10.1155/MI/2006/84829, LIF Upregulates Expression of NK-1R in NHBE Cells
2. PubMed ID: 33536035, Zhou L,Li C,Liu X,Zhang T.Effect of Irisin on LIF and integrin alphavbeta3 in rats of implantation failure.Reprod Biol Endocrinol.2021 Feb 3;19(1):18.doi:10.1186/s12958-021-00700-9.PMID:33536035;PMCID:PMC7856750.
3. PubMed ID: 17392578, Hu Cp, Feng Jt, Tang Yl, Zhu Jq, Lin Mj, Yu Me. Mediators Inflamm. 2006;2006(5):84829. Lif Upregulates Expression Of Nk-1R In Nhbe Cells.

Visit [bosterbio.com/anti-lif-picoband-trade-antibody-pb9036-boster.html](https://bosterbio.com/anti-lif-picoband-trade-antibody-pb9036-boster.html) to see all 10 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LIF Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.