

Anti-CNTF Antibody Picoband® PE Conjugated

Catalog Number: PB9067-PE

About Cntf

Ciliary neurotrophic factor is a protein that in humans is encoded by the CNTF gene. The protein encoded by this gene is a polypeptide hormone and nerve growth factor whose actions have mainly been studied in the nervous system where it promotes neurotransmitter synthesis and neurite outgrowth in certain neural populations including astrocytes. It is mapped to 11q12.1. In addition to neurotrophic effects on parasympathetic neurons, CNTF was shown to have activities on sympathetic and sensory neurons. CNTF appears not to be involved in motor neuron survival during development. Human CNTF has a trophic influence on degenerating striatal neurons as well as on critical nonstriatal regions. CNTF plays a role in a cardiac signal transduction pathway that regulates obesity-related left ventricular hypertrophy.

Overview

Product Name	Anti-CNTF Antibody Picoband® PE Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P20294

Technical Details

Immunogen	E. coli-derived rat CNTF recombinant protein (Position: A2-M200). Rat CNTF shares 84% and 95% amino acid (aa) sequences identity with human and mouse CNTF, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CNTF Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.