

## Anti-Hsp90 alpha/HSP90AA1 Antibody Picoband® FITC Conjugated

Catalog Number: PB9089-FITC

### About HSP90AA1

Heat shock protein HSP 90-alpha is a protein that in humans is encoded by the HSP90AA1 gene. The gene, HSP90AA1, encodes the human stress-inducible 90-kDa heat shock protein alpha (Hsp90A). Complemented by the constitutively expressed paralog Hsp90B which shares over 85% amino acid sequence identity, Hsp90A expression is initiated when a cell experiences proteotoxic stress. Once expressed Hsp90A dimers operate as molecular chaperones that bind and fold other proteins into their functional 3-dimensional structures. This molecular chaperoning ability of Hsp90A is driven by a cycle of structural rearrangements fueled by ATP hydrolysis. Current research on Hsp90A focuses in its role as a drug target due to its interaction with a large number of tumor promoting proteins and its role in cellular stress adaptation.

### Overview

Product Name	Anti-Hsp90 alpha/HSP90AA1 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P07900

### Technical Details

Immunogen	E.coli-derived human Hsp90 alpha recombinant protein (Position: P2-V365). Human Hsp90 alpha shares 99% amino acid (aa) sequence identity with both mouse and rat Hsp90 alpha.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 7 Publications Citing This Product

1. PubMed ID: 10.3390/ijms15010186, The Effect of 5-Adenylic Acid on Hepatic Proteome of Mice Radiated by 60Co gamma-ray
2. PubMed ID: 10.3892/ijo.2015.3039, Autophagy is involved in recombinant Newcastle disease virus (rL-RVG)-induced cell death of stomach adenocarcinoma cells in vitro
3. PubMed ID: 10.1016/j.cellbi.2009.01.008, Rat bone marrow derived mesenchymal progenitor cells support mouse ES cell growth and germ-like cell differentiation

Visit [bosterbio.com/anti-hsp90-alpha-picoband-trade-antibody-pb9089-boster.html](http://bosterbio.com/anti-hsp90-alpha-picoband-trade-antibody-pb9089-boster.html) to see all 7 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hsp90 alpha/HSP90AA1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.