

## Anti-c-Myc Antibody Picoband® Cy3 Conjugated

Catalog Number: PB9092-Cy3

### About MYC

C-Myc is an oncogene that functions both in the stimulation of cell proliferation and in apoptosis. C-Myc elicits its oncogenic activity by causing immortalization, and to a lesser extent the transformation of cells, in addition to several other mechanisms. The c-MYC proto-oncogene encodes a transcription factor that is critical for cell growth and proliferation. It is one of the genes frequently altered in cancer cells in which it exhibits constitutive activity. Downregulation of c-Myc is critical for 2-Methoxyestradiol (2ME2)-induced oxidative stress and apoptosis in AML cells. And its up-regulation is important for promoting lymphocyte cell division, and demonstrating that GFP-c-Myc expression is a marker of proliferating lymphocytes in vivo.

### Overview

Product Name	Anti-c-Myc Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P01106

### Technical Details

Immunogen	E.coli-derived human c-Myc recombinant protein (Position: E257-A439). Human c-Myc shares 91% amino acid (aa) sequences identity with both mouse and rat c-Myc.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 19 Publications Citing This Product

1. PubMed ID: 10.3748/wjg.14.5008, Positional and expressive alteration of prohibitin during the induced differentiation of human hepatocarcinoma SMMC-7721 cells
2. PubMed ID: 10.3892/or\_00000554, Let-7a microRNA functions as a potential tumor suppressor in human laryngeal cancer
3. PubMed ID: 24194897, Li M, Tian L, Wang L, Yao H, Zhang J, Lu J, Sun Y, Gao X, Xiao H, Liu M. Plos One. 2013 Oct 23;8(10):E77829. Doi: 10.1371/Journal.Pone.0077829. Ecollection 2013. Down-Regulation Of Mir-129-5P Inhibits Growth And Induces Apoptosis In Laryngeal Squa...

Visit [bosterbio.com/anti-c-myc-picoband-trade-antibody-pb9092-boster.html](http://bosterbio.com/anti-c-myc-picoband-trade-antibody-pb9092-boster.html) to see all 19 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-c-Myc Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.