

## Anti-Cyclin D3/CCND3 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9103-Fluoro488

### About CCND3

CCND3, also called Cyclin D3, is a protein that in humans is encoded by the CCND3 gene. It is mapped to 6p21.1. The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin forms a complex with and functions as a regulatory subunit of CDK4 or CDK6, whose activity is required for cell cycle G1/S transition. In addition, overexpression of CCND3 upregulated the translational activity in HeLa cells in a dose-dependent manner.

### Overview

Product Name	Anti-Cyclin D3/CCND3 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P30281

### Technical Details

Immunogen	E.coli-derived human Cyclin D3 recombinant protein (Position: Q136-L292). Human Cyclin D3 shares 94% and 93% amino acid (aa) sequences identity with mouse and rat Cyclin D3, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 10.1007/s10616-018-0235-3, Co-culture with TM4 cells enhances the proliferation and migration of rat adipose-derived mesenchymal stem cells with high stemness

Visit [bosterbio.com/anti-cyclin-d3-picoband-trade-antibody-pb9103-boster.html](http://bosterbio.com/anti-cyclin-d3-picoband-trade-antibody-pb9103-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cyclin D3/CCND3 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.