

## Anti-CIAS1/NALP3/NLRP3 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9114-Fluoro488

### About NLRP3

NLRP3 (NLR FAMILY, PYRIN DOMAIN-CONTAINING 3), also known as CIAS1, CRYOPYRIN, NALP3 or PYPAF1, is a protein that in humans is encoded by the NLRP3 gene. It is mapped to 1q44. The encoded protein may play a role in the regulation of inflammation and apoptosis. Mutation of the NALP3 nucleotide-binding domain reduced ATP binding, CASP1 activation, IL1B production, cell death, macromolecular complex formation, self-association, and association with ASC. Consistent with an essential role for Nlrp3 inflammasomes in antifungal immunity, it has been showed that Nlrp3-deficient mice are hypersusceptible to *C. albicans* infection. Activation of the NLRP3 inflammasome in response to virus or to RNA was dependent upon lysosomal maturation and reactive oxygen species production in human cells. The NLRP3 inflammasome senses obesity-associated danger signals and contributes to obesity-induced inflammation and insulin resistance.

### Overview

Product Name	Anti-CIAS1/NALP3/NLRP3 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96P20

### Technical Details

Immunogen	E.coli-derived human CIAS1 recombinant protein (Position: D21-K194). Human CIAS1 shares 77% amino acid (aa) sequence identity with mouse CIAS1.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

### 3 Publications Citing This Product

1. PubMed ID: 31156617, Wei Z,Wang J,Wang Y,Wang C,Liu X,Han Z,Fu Y,Yang Z.Effects of Neutrophil Extracellular Traps on Bovine Mammary Epithelial Cells in vitro.Front Immunol.2019 May 17;10:1003.doi:10.3389/fimmu.2019.01003.PMID:31156617;PMCID:PMC6533846.

2. PubMed ID: 24884473, Wang B, Yang G, Liang X, Zhu M, Du M. Bmc Complement Altern Med. 2014 May 20;14:162. Doi: 10.1186/1472-6882-14-162. Grape Seed Extract Prevents Skeletal Muscle Wasting In Interleukin 10 Knockout Mice.

3. PubMed ID: 24775094, Xue Y, Wang H, Du M, Zhu Mj. J Nutr Biochem. 2014 Jul;25(7):758-64. Doi: 10.1016/J.jnutbio.2014.03.009. Epub 2014 Apr 1. Maternal Obesity Induces Gut Inflammation And Impairs Gut Epithelial Barrier Function In Nonobese Diabetic Mice.

Visit [bosterbio.com/anti-cias1-nalp3-picoband-trade-antibody-pb9114-boster.html](http://bosterbio.com/anti-cias1-nalp3-picoband-trade-antibody-pb9114-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CIAS1/NALP3/NLRP3 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.