

Anti-AP2 alpha/TFAP2A Antibody Picoband® APC Conjugated

Catalog Number: PB9118-APC

About TFAP2A

TFAP2A, also known as AP-2alpha or BOFS, is a protein that in humans is encoded by the TFAP2A gene. It is mapped to 6p24.3. The protein encoded by this gene is a transcription factor that binds the consensus sequence 5'-GCCNNNGGC-3'. The encoded protein functions as either a homodimer or as a heterodimer with similar family members. This protein activates the transcription of some genes while inhibiting the transcription of others. TFAP2A acts as a sequence-specific DNA-binding transcription factor recognizing and binding to the specific DNA sequence and recruiting transcription machinery. This gene is expressed in neural crest cell lineages with the highest levels of expression corresponding to early neural crest cells, suggesting that TFAP2A plays a role in their differentiation and development.

Overview

| | |
|----------------------|--|
| Product Name | Anti-AP2 alpha/TFAP2A Antibody Picoband® APC Conjugated |
| Reactive Species | Human, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P05549 |

Technical Details

| | |
|------------------|---|
| Immunogen | E.coli-derived human AP2 alpha recombinant protein (Position: M1-G166). Human AP2 alpha shares 98% amino acid (aa) sequence identity with both mouse and rat AP2 alpha. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AP2 alpha/TFAP2A Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.