

Anti-AP2 alpha/TFAP2A Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9118-Fluoro550

About TFAP2A

TFAP2A, also known as AP-2alpha or BOFS, is a protein that in humans is encoded by the TFAP2A gene. It is mapped to 6p24.3. The protein encoded by this gene is a transcription factor that binds the consensus sequence 5'-GCCNNNGGC-3'. The encoded protein functions as either a homodimer or as a heterodimer with similar family members. This protein activates the transcription of some genes while inhibiting the transcription of others. TFAP2A acts as a sequence-specific DNA-binding transcription factor recognizing and binding to the specific DNA sequence and recruiting transcription machinery. This gene is expressed in neural crest cell lineages with the highest levels of expression corresponding to early neural crest cells, suggesting that TFAP2A plays a role in their differentiation and development.

Overview

Product Name	Anti-AP2 alpha/TFAP2A Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P05549

Technical Details

Immunogen	E.coli-derived human AP2 alpha recombinant protein (Position: M1-G166). Human AP2 alpha shares 98% amino acid (aa) sequence identity with both mouse and rat AP2 alpha.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AP2 alpha/TFAP2A Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.