

Anti-AP2 alpha/TFAP2A Antibody Picoband® HRP Conjugated

Catalog Number: PB9118-HRP

About TFAP2A

TFAP2A, also known as AP-2alpha or BOFS, is a protein that in humans is encoded by the TFAP2A gene. It is mapped to 6p24.3. The protein encoded by this gene is a transcription factor that binds the consensus sequence 5'-GCCNNNGGC-3'. The encoded protein functions as either a homodimer or as a heterodimer with similar family members. This protein activates the transcription of some genes while inhibiting the transcription of others. TFAP2A acts as a sequence-specific DNA-binding transcription factor recognizing and binding to the specific DNA sequence and recruiting transcription machinery. This gene is expressed in neural crest cell lineages with the highest levels of expression corresponding to early neural crest cells, suggesting that TFAP2A plays a role in their differentiation and development.

Overview

Product Name	Anti-AP2 alpha/TFAP2A Antibody Picoband® HRP Conjugated
Reactive Species	Human, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P05549

Technical Details

Immunogen	E.coli-derived human AP2 alpha recombinant protein (Position: M1-G166). Human AP2 alpha shares 98% amino acid (aa) sequence identity with both mouse and rat AP2 alpha.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users.

Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.
ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AP2 alpha/TFAP2A Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.