

Anti-Decorin/DCN Antibody Picoband® HRP Conjugated

Catalog Number: PB9174-HRP

About DCN

Decorin is a protein that in humans is encoded by the DCN gene. This gene is mapped to 12q21.3. It belongs to the small leucine-rich proteoglycan (SLRP) family and consists of a protein core containing leucine repeats with a glycosaminoglycan (GAG) chain consisting of either chondroitin sulfate (CS) or dermatan sulfate (DS). Decorin is a small cellular or pericellular matrix proteoglycan and is closely related in structure to biglycan protein. This protein is a component of connective tissue, binds to type I collagen fibrils, and plays a role in matrix assembly. And it also may play a role in epithelial/mesenchymal interactions during organ development and shaping. Decorin has been shown to have anti-tumorigenic properties in an experimental murine tumor model and is capable of suppressing the growth of various tumor cell lines.

Overview

Product Name	Anti-Decorin/DCN Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P07585

Technical Details

Immunogen	E.coli-derived human Decorin recombinant protein (Position: D31-K359). Human Decorin shares 80% and 77% amino acid (aa) sequences identity with mouse and rat Decorin, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

3 Publications Citing This Product

1. PubMed ID: 27713567, An Animal Model of Modic Changes by Embedding Autogenous Nucleus Pulposus inside Subchondral Bone of Lumbar Vertebrae
2. PubMed ID: 24987434, Huang Gf, Zou J, Shi J, Zhang Dy, Pen Hf, Zhang Q, Gao Y, Wang By. Evid Based Complement Alternat Med. 2014;2014:731395. Doi: 10.1155/2014/731395. Epub 2014 May 27. The Effect Of Electroacupuncture On The Extracellular Matrix Synthesis And Degrada...
3. PubMed ID: 28290552, Lai J, Chen F, Chen J, Ruan G, He M, Chen C, Tang J, Wang DW. Sci Rep. 2017 Mar 14;7:44473. doi: 10.1038/srep44473. Overexpression of decorin promoted angiogenesis in diabetic cardiomyopathy via IGF1R-AKT-VEGF signaling

Visit bosterbio.com/anti-decorin-picoband-trade-antibody-pb9174-boster.html to see all 3 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Decorin/DCN Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.