

## Anti-C5/C5a Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9187-Fluoro488

### About C5

C5a is a protein fragment released from complement component C5. This gene is mapped to 9q33.2. The protein encoded by this gene is the fifth component of complement, which plays an important role in inflammatory and cell killing processes. This protein is comprised of alpha and beta polypeptide chains that are linked by a disulfide bridge. An activation peptide, C5a, which is an anaphylatoxin that possesses potent spasmogenic and chemotactic activity, is derived from the alpha polypeptide via cleavage with a convertase. C5a plays a key role in increasing migration and adherence of neutrophils and monocytes to vessel walls. Mutations in this gene cause complement component 5 deficiency, a disease where patients show a propensity for severe recurrent infections. Defects in this gene have also been linked to susceptibility to liver fibrosis and to rheumatoid arthritis.

### Overview

Product Name	Anti-C5/C5a Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P08650

### Technical Details

Immunogen	E.coli-derived rat C5a recombinant protein (Position: D1-R77). Rat C5a shares 62% and 82% amino acid (aa) sequences identity with human and mouse C5a, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm

	Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 29956782, C5a aggravates dysfunction of the articular cartilage and synovial fluid in rats with knee joint immobilization

Visit [bosterbio.com/anti-c5a-picoband-trade-antibody-pb9187-boster.html](https://bosterbio.com/anti-c5a-picoband-trade-antibody-pb9187-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C5/C5a Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.