

Anti-HOXA10 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9217-Fluoro488

About HOXA10

Homeobox protein Hox-A10 is a protein that in humans is encoded by the HOXA10 gene. HOXA10 is expressed in the adult human endometrium, and the expression of HOXA10 dramatically increased during the midsecretory phase of the menstrual cycle, corresponding to the time of implantation and increase in circulating progesterone. HOXA10 may have an important function in regulating endometrial development during the menstrual cycle and in establishing conditions necessary for implantation in the human. HOXA10 is also a master regulator of postnatal hematopoietic development whose tight regulation is required for normal development of erythroid and megakaryocytic lineages.

Overview

Product Name	Anti-HOXA10 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P31260

Technical Details

Immunogen	E.coli-derived human HOXA10 recombinant protein (Position: Y57-R340). Human HOXA10 shares 87% amino acid (aa) sequence identity with mouse HOXA10.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HOXA10 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.