

Anti-IKK beta/IKBKB Antibody Picoband® HRP Conjugated

Catalog Number: PB9221-HRP

About IKBKB

IKBKB (Inhibitor of Kappa Light Chain Gene Enhancer in B Cells, Kinase of, Beta), also known as IKKB or NFKB1KB, is a protein that in humans is encoded by the IKBKB gene. This gene is mapped to chromosome 8p12-p11 by FISH. It is found that mutations in IKK2 had a more pronounced effect upon NFkB activation than did comparable mutations in IKK1. The activity of various antiinflammatory agents on the IKK complex has been tested, it has been found that aspirin and sodium salicylate specifically inhibit IKK-beta activity in vitro and in vivo by binding to IKK-beta to reduce ATP binding.

Overview

Product Name	Anti-IKK beta/IKBKB Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O14920

Technical Details

Immunogen	E.coli-derived human IKK beta recombinant protein (Position: E398-S756). Human IKK beta shares 91% and 90% amino acid (aa) sequences identity with mouse and rat IKK beta, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IKK beta/IKBKB Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.