

## Anti-MUM1/IRF4 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9222-Fluoro550

### About IRF4

Interferon regulatory factor 4 (IRF4), also known as MUM1, is a protein that in humans is encoded by the IRF4 gene. It is located on 6p25.3. IRF4 is a transcription factor, and it is essential for the development of T helper-2 (Th2) cells, IL17 -producing Th17 cells, and IL9 -producing Th9 cells. In melanocytic cells, the IRF4 gene may be regulated by MITF. IRF4 is a transcription factor that has been implicated in acute leukemia. This gene is strongly associated with pigmentation, sensitivity of skin to sun exposure, freckles, blue eyes, and brown hair color. What's more, IRF4 inhibition is toxic to myeloma cell lines, regardless of transforming oncogenic mechanism.

### Overview

Product Name	Anti-MUM1/IRF4 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q15306

### Technical Details

Immunogen	E.coli-derived human MUM1 recombinant protein (Position: E272-E451). Human MUM1 shares 92% amino acid (aa) sequence identity with mouse MUM1.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

---

## 1 Publications Citing This Product

---

1. PubMed ID: 10.3892/or.2016.4562, TIPE2, a negative regulator of TLR signaling, regulates p27 through IRF4-induced signaling

Visit [bosterbio.com/anti-mum1-picoband-trade-antibody-pb9222-boster.html](http://bosterbio.com/anti-mum1-picoband-trade-antibody-pb9222-boster.html) to see all 1 publications.

---

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MUM1/IRF4 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.