

## Anti-IRS1 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9223-Fluoro647

### About IRS1

Insulin receptor substrate 1 (IRS-1) is a signalling adapter protein that in humans is encoded by the IRS-1 gene. It is mapped to 2q36.3. This gene exhibited no intrinsic enzyme activity, and it can serve as a docking protein involved in binding and activating other signal transduction molecules after being phosphorylated on tyrosine by insulin receptor kinase. IRS1 plays a key role in transmitting signals from the insulin and insulin-like growth factor-1 (IGF-1) receptors to intracellular pathways PI3K/Akt and Erk MAP kinase pathways. IRS1 also has important biological function for both metabolic and mitogenic (growth promoting) pathways. In addition to those, IRS1 is a key regulator of PI3K within malignant cells.

### Overview

Product Name	Anti-IRS1 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P35568

### Technical Details

Immunogen	E.coli-derived human IRS1 recombinant protein (Position: S1041-Q1242). Human IRS1 shares 78% and 80% amino acid (aa) sequences identity with mouse and rat IRS1, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

1. PubMed ID: PMID:26550236, Vibration exercise decreases insulin resistance and modulates the insulin signaling pathway in a type 2 diabetic rat model

2. PubMed ID: 25536154, Mi L, Chen Y, Zheng X, Li Y, Zhang Q, Mo D, Yang G. J Cell Biochem. 2015 Jul;116(7):1195-204. Doi: 10.1002/Jcb.25065. Microrna-139-5P Suppresses 3T3-L1 Preadipocyte Differentiation Through Notch And Irs1/Pi3K/Akt Insulin Signaling Pathways.

Visit [bosterbio.com/anti-irs1-picoband-trade-antibody-pb9223-boster.html](http://bosterbio.com/anti-irs1-picoband-trade-antibody-pb9223-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IRS1 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.