

Anti-Kv4.3/KCND3 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9229-Fluoro488

About KCND3

Potassium voltage-gated channel subfamily D member 3, also known as Kv4.3, is a protein that in humans is encoded by the KCND3 gene. KCND3 is a member of the potassium channel, voltage-gated, shal-related subfamily. Voltage-gated potassium (Kv) channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. It is mapped to 1p13.2. KCND3 is important in membrane repolarization in excitable cells. It contributes to the cardiac transient outward potassium current (Ito1), the main contributing current to the repolarizing phase 1 of the cardiac action potential.

Overview

Product Name	Anti-Kv4.3/KCND3 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9UK17

Technical Details

Immunogen	E.coli-derived human Kv4.3 recombinant protein (Position: M1-H177). Human Kv4.3 shares 100% and 99% amino acid (aa) sequences identity with mouse and rat Kv4.3, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Kv4.3/KCND3 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.